

Hydrophilic Lens :

- Vistazo is made of a hydrophilic acrylic material, with proven reputation and many years of clinical experience.
- Vistazo is characterized by excellent biocompatibility and mechanical quality.
- Vistazo material incorporates a natural yellow, violet filtering, chromophore for better protection of the retina.
- Ensures excellent stability and centration due to the unique C-loop mechanical design of the haptics.
- 360° square edge in order to minimize PCO.



MODEL NO.-CL160125CR6A

C-Loop

KEY SPECIFICATIONS

Optic Diameter	6.00 mm
Overall Diameter	12.50 mm
Haptic Design	C-Loop



MODEL NO.-LL160130CR5A

L-Loop

KEY SPECIFICATIONS

Optic Diameter	6.00 mm
Overall Diameter	13.00 mm
Haptic Design	L-Loop



MODEL NO.-FL160125CR6A

Flex Loop

KEY SPECIFICATIONS

Optic Diameter	6.00 mm
Overall Diameter	12.50 mm
Haptic Design	Flex Loop



Parameters	Ultrasound Biometry	Optical Biometry
A constant:SRK T	118	118.4
A constant:SRK II	118.2	118.65
Holladay I (Sf)	1.12	1.35
Hoffer Q	pACD-4.91	pACD-5.14
Haigis	a0-0.673	a0-0.912
	a1-0.400	a1-0.400
	a2-0.100	a2-0.100
Barrett LF		1.57

Optic Design - Aspheric/Bi-Convex, o Edge Design - 360° Continuous Square Edge, Haptic Angle - 0° / 5°, Diopter Range : +5D to +30D(15D to 26D in 0.5 D increments) *Available In Clear, Extended & Natural Yellow
**Other Diopters Available On Request

Available in Preloaded System

MODEL NO.-QL160110CR5A /
QL160110EY6A / QL160110NY5A

Monofocal

KEY SPECIFICATIONS

Optic Diameter

6.00 mm

Overall Diameter

11.00 mm

Haptic Design

Quadra Loop



MODEL NO.-QL160110NY5Tx

Toric

KEY SPECIFICATIONS

Optic Diameter

6.00 mm

Overall Diameter

11.00 mm

Haptic Design

Quadra Loop



MODEL NO.-QL160110NY5M

Trifocal

KEY SPECIFICATIONS

Optic Diameter

6.00 mm

Overall Diameter

11.00 mm

Haptic Design

Quadra Loop



**ENGINEER
THE POWER OF
CHANGE**

OPHTHALMIC
Solution, Inc.
Manufacturer of Intraocular Lenses

Address : 820 Corridor Park Blvd. Knoxville, TN 37932

Phone : (P) +1-865-671-4300 (F) +1-865-671-4302

Email : info@osilens.com

Website : www.opthalmicsolutioninc.com